

Illustrated Dyes	Concentration (g/l)	pH	Light Fastness
Specialty Gold N (F.A.O.)	15.0	4.5	8
Specialty Gold S	2.0	6.0	5
Specialty Fast Gold LN	2.0	5.5	7
Specialty Gold EAN	2.0	6.5	4
Specialty Golden Orange RLW	2.0	7.0	7
Specialty Orange 3A	2.0	6.0	6
Specialty Orange 2B	2.0	6.0	6
Specialty Yellow D	2.0	6.0	5
Specialty Yellow 4A	2.0	6.0	6
Specialty Yellow 3GL	3.0	6.0	8
Specialty Bronze 2LN	2.0	5.5	7
Specialty Copper BF	2.0	7.0	3
Specialty Brown BL	2.0	6.5	5
Specialty Brown GL	2.0	6.5	5
Specialty Bordeaux R	2.0	6.0	3
Specialty Deep Red L	5.0	5.5	6
Specialty Neon Pink	2.0	5.5	6
Specialty Violet 3D	2.0	6.0	3
Specialty Violet DS	2.0	6.0	3
Specialty Violet MRB	2.0	5.5	4
Specialty Fast Blue 2LWN	2.0	5.5	7
Specialty Fast Blue GN	4.0	5.7	8
Specialty Blue 4A	2.0	6.0	4
Specialty Turquoise PLWN	5.0	5.5	8
Specialty Blue A	2.0	6.0	8
Specialty Teal	4.0	5.5	7
Specialty Green AEN	2.0	6.0	7
Specialty Green SCG	4.0	5.5	6
Specialty Olive Drab	2.0	7.5	3
Specialty Grey NLN	2.0	5.5	6
Specialty Grey BL	2.0	5.5	6
Specialty Black BK Excel	10.0	5.5	7
Specialty Deep Black HL3 Paste	10.0	5.5	8